

Phospho-TrkB (Tyr705) Rabbit mAb


ABZ4351

Product Overview

Name	Phospho-TrkB (Tyr705) Rabbit mAb
Accession(Primary)	Q16620
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 140,90 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16620)

Gene Name	NTRK2
About Gene Symbol	



Western blot analysis of Phospho-TrkB (Tyr705) in HeLa, A549, C6 lysates using Phospho-TrkB (Tyr705) antibody.