

Phospho-TrkA/B (Tyr490/Tyr516) Rabbit mAb

ABZ4352


Product Overview

Name	Phospho-TrkA/B (Tyr490/Tyr516) Rabbit mAb
Accession(Primary)	Q16620
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 92 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q16620)

Gene Name NTRK2

About Gene Symbol



Western blot analysis of Phospho-TrkA (tyr490)/TrkB (Tyr516) in K562, C6, 3T3 lysates using Phospho-TrkA/B (Tyr490/Tyr516) antibody.

