

Phospho-CDC37 (Ser13) Rabbit mAb


ABZ4353

Product Overview


Name	Phospho-CDC37 (Ser13) Rabbit mAb
Accession(Primary)	Q16543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q16543)

Gene Name	CDC37
About Gene Symbol	CDC37 (Cell Division Cycle 37, HSP90 Cochaperone) is a molecular chaperone involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellu



Western blot analysis of Phospho-Cdc37 (Ser13) in mouse liver lysates using Phospho-CDC37 (Ser13) antibody.



Western blot analysis of Cdc37 (Phospho-Ser13) in HeLa, rat Brain, Ramos lysates using Cdc37 (Phospho-Ser13) antibody.