

Phospho-p38 MAPK (Thr180) Rabbit mAb


ABZ4354

Product Overview

Name	Phospho-p38 MAPK (Thr180) Rabbit mAb
Accession(Primary)	Q16539
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; FC-1:200-1:500

Target information (Q16539)

Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha



Western blot analysis of Phospho-P38 (Thr180) in 3T3, HeLa lysates using Phospho-p38 (Thr180) antibody.