

Phospho-NRF2 (Ser40) Rabbit mAb


ABZ4355

Product Overview


Name	Phospho-NRF2 (Ser40) Rabbit mAb
Accession(Primary)	Q16236
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:500; ICC/IF-1:100-1:200; FC-1:50-1:100

Target information (Q16236)

Gene Name	NFE2L2
About Gene Symbol	



Western blot analysis of Phospho-Nrf2 (Ser40) in K562, HeLa lysates using Phospho-Nrf2 (Ser40) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Nrf2 (Phospho- S40) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.