

Phospho-Smad2 (Ser250) Rabbit mAb


ABZ4359

Product Overview


Name	Phospho-Smad2 (Ser250) Rabbit mAb
Accession(Primary)	Q15796
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 58,62 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:20-1:50

Target information (Q15796)

Gene Name	SMAD2
About Gene Symbol	SMAD2 (SMAD Family Member 2) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Smad2 (Ser250) in mouse brain, mouse lung lysates using Phospho-Smad2 (Ser250) antibody.



Western blot analysis of Smad2 (Phospho-Ser250) in HeLa, CHO-K1, C6 lysates using Smad2 (Phospho-Ser250) antibody.