

Phospho-Smad2/3 (Thr8) Rabbit pAb


ABZ4360

Product Overview


Name	Phospho-Smad2/3 (Thr8) Rabbit pAb
Accession(Primary)	Q15796
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 52 kDa; Observed MW: 48-52□60 kDa
Recommended Dilution	WB-1:500-1:2000; IF/ICC1:100-1:500

Target information (Q15796)

Gene Name	SMAD2/SMAD3
About Gene Symbol	SMAD2 (SMAD Family Member 2) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Smad2/3 (Thr8) in various lysates using Phospho-Smad2/3 (Thr8) antibody.



Western blot analysis of Phospho-Smad2/3 (Thr8) in VEC lysates using Phospho-Smad2/3 (Thr8) antibody.