


Phospho-ErbB 4 (Tyr1162) Rabbit mAb

ABZ4362


Product Overview

Name	Phospho-ErbB 4 (Tyr1162) Rabbit mAb
Accession(Primary)	Q15303
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 147 kDa; Observed MW: 180 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; IP-1:20

Target information (Q15303)

Gene Name ERBB4**About Gene Symbol**

Western blot analysis of Phospho-ErbB 4 (Tyr1162) in HeLa, A549 lysates using Phospho-ErbB 4 (Tyr1162) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using ErbB 4 (Phospho-Tyr1162) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.