

Phospho-KIF22 (Ser427) Rabbit mAb


ABZ4363

Product Overview

Name	Phospho-KIF22 (Ser427) Rabbit mAb
Accession(Primary)	Q14807
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q14807)

Gene Name	KIF22
About Gene Symbol	KIF22 (Kinesin Family Member 22) is a kinesin motor protein involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-KIF22 (Ser427) in A549, HL-60, U2OS lysates using Phospho-KIF22 (Ser427) antibody.