

Phospho-PTK2B (Tyr402) Rabbit pAb

ABZ4364

Product Overview

Name	Phospho-PTK2B (Tyr402) Rabbit pAb
Accession(Primary)	Q14289
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:300; ELISA-1:20000; IF-1:50-1:200

Target information (Q14289)

Gene Name	PTK2B
About Gene Symbol	

