

Phospho-4E BP1 (Thr70) Rabbit pAb


ABZ4366

Product Overview

Name	Phospho-4E BP1 (Thr70) Rabbit pAb
Accession(Primary)	Q13541
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q13541)

Gene Name	EIF4EBP1
About Gene Symbol	EIF4EBP1 (Eukaryotic Translation Initiation Factor 4E Binding Protein 1) is a protein kinase and signaling regulator involved in mTOR/PI3K/Akt signaling and translational regulation. It is an important research target in cancer, cell signaling, and therap



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma tissue using Phospho-4E BP1 (Thr70) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.