

Phospho-AMPK alpha 1 (Ser496) Rabbit mAb


ABZ4368

Product Overview


Name	Phospho-AMPK alpha 1 (Ser496) Rabbit mAb
Accession(Primary)	Q13131
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:100-1:200

Target information (Q13131)

Gene Name	PRKAA1
About Gene Symbol	PRKAA1 (Protein Kinase AMP-Activated Catalytic Subunit Alpha 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is



Western blot analysis of AMPK alpha 1 (Phospho-Ser496) in Ramos lysates using Phospho-AMPK alpha 1 (Ser496) antibody.



Immunocytochemistry analysis of AMPK alpha 1 (Phospho-Ser496) (green) in HeLa using AMPK alpha 1 (Phospho-Ser496) antibody, and DAPI (blue)