

## Phospho-FOXO1/3 (Ser322/S325) Rabbit pAb


ABZ4370

### Product Overview

|                             |                                                                  |
|-----------------------------|------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-FOXO1/3 (Ser322/S325) Rabbit pAb                         |
| <b>Accession(Primary)</b>   | Q12778                                                           |
| <b>Host Species</b>         | Rabbit                                                           |
| <b>Clonality</b>            | Polyclonal                                                       |
| <b>Isotype</b>              | IgG                                                              |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                |
| <b>Calculated MW</b>        | Calculated MW: 70 kDa; Observed MW: 97 kDa                       |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:5000 |

### Target information (Q12778)

|                          |                                                                                                                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FOXO1/FOXO3                                                                                                                                                                                                                                     |
| <b>About Gene Symbol</b> | FOXO1 (Forkhead Box O1) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Phospho-FOXO1/3 (Ser322/S325) in 3T3 lysates using Phospho-FOXO1/3 (Ser322/S325) antibody.