

Phospho-MEK1 (Thr292) Rabbit mAb


ABZ4377

Product Overview

Name	Phospho-MEK1 (Thr292) Rabbit mAb
Accession(Primary)	Q02750
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q02750)

Gene Name	MAP2K1
About Gene Symbol	MAP2K1 (Mitogen-Activated Protein Kinase Kinase 1) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-MEK1 (Thr292) in C6, 3T3 lysates using Phospho-MEK1 (Thr292) antibody.