

Phospho-PKC epsilon (Ser729) Rabbit pAb

ABZ4378

Product Overview

Name	Phospho-PKC epsilon (Ser729) Rabbit pAb
Accession(Primary)	Q02156
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q02156)

Gene Name

PRKCE

About Gene Symbol

 PKC epsilon --
 (pSer729)


-- 117

-- 85

-- 48


-- 34

-- 26

-- 19

(kD)

Western blot analysis of Phospho-PKC epsilon (Ser729) in HeLa lysates using Phospho-PKC epsilon (Ser729) antibody.



Immunohistochemistry analysis of paraffin-embedded Human Liver Cancer using Phospho-PKC epsilon (Ser729) antibody. High-pressure and temperature Tris-EDTA pH 9.0 was used for antigen retrieval.