

Phospho-MDM2 (Ser166) Rabbit mAb

ABZ4379


Product Overview


Name	Phospho-MDM2 (Ser166) Rabbit mAb
Accession(Primary)	Q00987
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (Q00987)

Gene Name MDM2

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Phospho-MDM2 (Ser166) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.