

Phospho-CDK5 (Tyr15) Rabbit pAb


ABZ4380

Product Overview

Name	Phospho-CDK5 (Tyr15) Rabbit pAb
Accession(Primary)	Q00535
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q00535)

Gene Name	CDK5
About Gene Symbol	CDK5 (Cell division protein kinase 5) is a protein kinase and signaling regulator involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-CDK5 (Tyr15) in COS7 lysates using Phospho-CDK5 (Tyr15) antibody.