

Phospho-Smad3 (Ser425) (7H6) Mouse mAb


ABZ4381

Product Overview


Name	Phospho-Smad3 (Ser425) (7H6) Mouse mAb
Accession(Primary)	P84022
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100

Target information (P84022)

Gene Name	SMAD3
About Gene Symbol	SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Skin Tissue using Phospho-Smad3 (Ser425) (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Phospho-Smad3 (Ser425) (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.