

Phospho-Smad3 (Ser423/425) Rabbit mAb


ABZ4383

Product Overview

Name	Phospho-Smad3 (Ser423/425) Rabbit mAb
Accession(Primary)	P84022
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 48 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-1:50-1:200; ICC/IF-1:100-1:200

Target information (P84022)

Gene Name	SMAD3
About Gene Symbol	SMAD3 (SMAD Family Member 3) is a transcription factor involved in TGF- β signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-Smad3
(Ser423/425) in Jurkat, C6, CHO-K1, HeLa lysates
using Phospho-Smad3 (Ser423/425) antibody.