

Phospho-Histone H3 (Thr3) Rabbit mAb


ABZ4385

Product Overview


Name	Phospho-Histone H3 (Thr3) Rabbit mAb
Accession(Primary)	P68431
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:100

Target information (P68431)

Gene Name	H3C1
About Gene Symbol	H3C1 (H3 Clustered Histone 1) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to el



Western blot analysis of Phospho-Histone H3 (Thr3) in C6, 3T3, HeLa lysates using Phospho-Histone H3 (Thr3) antibody.



Immunocytochemistry analysis of Phospho-Histone H3 (Thr3)(green) in HeLa using Phospho-Histone H3 (Thr3) antibody, and DAPI(blue)