

Phospho-Rac1 (Ser71) Rabbit pAb

ABZ4388

Product Overview

Name	Phospho-Rac1 (Ser71) Rabbit pAb
Accession(Primary)	P63000
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P63000)

Gene Name RAC1

About Gene Symbol



