

Phospho-LIM Kinase 1 (Thr508) Rabbit pAb


ABZ4389

Product Overview


Name	Phospho-LIM Kinase 1 (Thr508) Rabbit pAb
Accession(Primary)	P53667
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (P53667)


Gene Name	LIMK1
About Gene Symbol	



Western blot analysis of Phospho-LIM Kinase 1 (Thr508) in C2C12, HeLa, 293 lysates using Phospho-LIM Kinase 1 (Thr508) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using LIMK1(Phospho-Thr508) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of Phospho-LIM Kinase 1 (Thr508) in HeLa cells using LIMK1(Phospho-Thr508) antibody (red).