

Phospho-RSK2 (Ser227) Rabbit mAb

ABZ4391


Product Overview

Name	Phospho-RSK2 (Ser227) Rabbit mAb
Accession(Primary)	P51812
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:1000


Target information (P51812)

Gene Name RPS6KA3

About Gene Symbol



Western blot analysis of Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) in C6, C6+UV lysates using Phospho-RSK2 (Ser227) antibody.



Western blot analysis of Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) in 3T3, HeLa lysates using Phospho-Rsk 2/MAPKAP Kinase 1b (Ser227) antibody.