

Phospho-GSK3(alpha + beta)(Tyr279/Tyr216) Rabbit mAb


ABZ4393

Product Overview

Name	Phospho-GSK3(alpha + beta)(Tyr279/Tyr216) Rabbit mAb
Accession(Primary)	P49840
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 51 kDa; Observed MW: 47-51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; IP-1:20-1:50

Target information (P49840)

Gene Name	GSK3A/GSK3B
About Gene Symbol	GSK3A (Glycogen Synthase Kinase 3 Alpha) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-GSK3
(alpha/beta) (Tyr216/Tyr279) in C6, 3T3 lysates
using Phospho-GSK3 (Tyr216/Tyr279) antibody.