

Phospho-MCM2 (Ser41) Rabbit mAb

ABZ4394


Product Overview

Name	Phospho-MCM2 (Ser41) Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20


Target information (P49736)

Gene Name MCM2

About Gene Symbol



Western blot analysis of Phospho-MCM2 (Ser41) in Jurkat, C6, CHO-K1, HeLa lysates using Phospho-MCM2 (Ser41) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM2 (Phospho-S41) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.