


Phospho-MCM2 (Ser108) Rabbit mAb

ABZ4395


Product Overview

Name	Phospho-MCM2 (Ser108) Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:20-1:50; IP-1:20


Target information (P49736)

Gene Name MCM2**About Gene Symbol**

Western blot analysis of Phospho-MCM2 (Ser108) in rat spleen lysates using Phospho-MCM2 (Ser108) antibody.



Western blot analysis of Phospho-MCM2 (Ser108) in K562, 3T3, HeLa lysates using Phospho-MCM2 (Ser108) antibody.



Immunocytochemistry analysis of MCM2 (Phospho-Ser108)(green) in HeLa using MCM2 (Phospho-Ser108) antibody, and DAPI(blue)