

## Phospho-MCM2 (Ser27) Rabbit mAb


ABZ4396

### Product Overview


|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-MCM2 (Ser27) Rabbit mAb                                           |
| <b>Accession(Primary)</b>   | P49736                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                         |
| <b>Calculated MW</b>        | Calculated MW: 102 kDa; Observed MW: 125 kDa                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200; IP-1:10-1:100 |

### Target information (P49736)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | MCM2 |
| <b>About Gene Symbol</b> |      |



Immunocytochemistry analysis of Phospho-MCM2 (Ser27) (green) in HeLa using Phospho-MCM2 (Ser27) antibody, and DAPI (blue).



Western blot analysis of MCM2 (Phospho-Ser27) in HeLa, A549, HL-60, U2OS lysates using MCM2 (Phospho-Ser27) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM2 (Phospho-Ser27) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.