

Phospho-MCM2 (Ser40) Rabbit mAb

ABZ4397


Product Overview


Name	Phospho-MCM2 (Ser40) Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P49736)

Gene Name MCM2

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM2 (Phospho-Ser40) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.