

Phospho-YAP1 (Ser127) Rabbit pAb


ABZ4398

Product Overview


Name	Phospho-YAP1 (Ser127) Rabbit pAb
Accession(Primary)	P46937
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (P46937)


Gene Name	YAP1
About Gene Symbol	YAP1 (Yes1 Associated Transcriptional Regulator) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-YAP (Ser127) in COS7 lysates using Phospho-YAP (Ser127) antibody.



Immunohistochemistry analysis of paraffin-embedded Human heart using YAP (Phospho-Ser127) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of Phospho-YAP (Ser127) in HeLa cells using YAP (Phospho-Ser127) antibody (green) and beta-tubulin antibody (red).