

Phospho-YAP1 (Ser127) Rabbit mAb


ABZ4399

Product Overview

Name	Phospho-YAP1 (Ser127) Rabbit mAb
Accession(Primary)	P46937
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 70-75 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200

Target information (P46937)

Gene Name	YAP1
About Gene Symbol	YAP1 (Yes1 Associated Transcriptional Regulator) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-YAP1 (Ser127) in HeLa, HeLa+CalyculinA lysates using Phospho-YAP1 (Ser127) antibody.