

Phospho-MEK4 (Ser80) Rabbit mAb

ABZ4401


Product Overview


Name	Phospho-MEK4 (Ser80) Rabbit mAb
Accession(Primary)	P45985
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100

Target information (P45985)

Gene Name MAP2K4

About Gene Symbol





Western blot analysis of Phospho-MEK4/MKK4 (Ser80) in A549, MCF-7, U251 lysates using Phospho-MEK4/MKK4 (Ser80) antibody.