

## Phospho-JNK1 (Thr183/Tyr185) Rabbit mAb


ABZ4403

### Product Overview

|                             |                                               |
|-----------------------------|-----------------------------------------------|
| <b>Name</b>                 | Phospho-JNK1 (Thr183/Tyr185) Rabbit mAb       |
| <b>Accession(Primary)</b>   | P45983                                        |
| <b>Host Species</b>         | Rabbit                                        |
| <b>Clonality</b>            | Monoclonal                                    |
| <b>Isotype</b>              | IgG                                           |
| <b>Species Reactivity</b>   | Human                                         |
| <b>Calculated MW</b>        | Calculated MW: 48 kDa; Observed MW: 46,54 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20-1:50                 |

### Target information (P45983)

|                          |                                                                                                                                                                                                                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MAPK8                                                                                                                                                                                                                                                                                                   |
| <b>About Gene Symbol</b> | MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Phospho-JNK1  
(Thr183/Tyr185) in HeLa, A549, C6 lysates using  
Phospho-JNK1 (Thr183/Tyr185) antibody.