

Phospho-JNK1/2/3 (Thr183/Tyr185) Rabbit pAb


ABZ4404

Product Overview

Name	Phospho-JNK1/2/3 (Thr183/Tyr185) Rabbit pAb
Accession(Primary)	P45983
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Chicken, Pig
Calculated MW	Calculated MW: 48 kDa; Observed MW: 46,54 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:200; IHC-1:50-1:100

Target information (P45983)

Gene Name	MAPK8
About Gene Symbol	MAPK8 (Mitogen-Activated Protein Kinase 8) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-JNK
(Thr183/Tyr185) in various lysates using Phospho-
p-JNK1/2/3 (T183/Y185) JNK (Thr183/Tyr185) antibody.