

Phospho-Syk (Tyr348) Rabbit pAb


ABZ4405

Product Overview


Name	Phospho-Syk (Tyr348) Rabbit pAb
Accession(Primary)	P43405
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P43405)

Gene Name	SYK
About Gene Symbol	



Western blot analysis of Phospho-Syk (Tyr348) in 3T3 lysates using Phospho-Syk (Tyr348) antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Phospho-Syk (Tyr348) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.