

Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845)


ABZ4408

Product Overview


Name	Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845) Rabbit mAb
Accession(Primary)	P42261
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 102 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P42261)

Gene Name	GRIA1
About Gene Symbol	



Western blot analysis of Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) in U87-MG lysates using Phospho-Glutamate Receptor 1 (AMPA Subtype) (Ser845) antibody.



Western blot analysis of Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) in C6 lysates using Phospho-Glutamate Receptor 1 (AMPA subtype) (Ser845) antibody.