

Phospho-STAT1 (Ser727) Rabbit mAb


ABZ4410

Product Overview


Name	Phospho-STAT1 (Ser727) Rabbit mAb
Accession(Primary)	P42224
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:200-1:1000

Target information (P42224)

Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell



Western blot analysis of Phospho-STAT1 (Ser727) in HeLa, HeLa+UV, C6, C6+UV lysates using Phospho-STAT1 (Ser727) antibody.



Western blot analysis of Phospho-STAT1 (Ser727) in 3T3 lysates using Phospho-STAT1 (Ser727) antibody.