

Phospho-STAT1 (Tyr701) Rabbit pAb


ABZ4411

Product Overview

Name	Phospho-STAT1 (Tyr701) Rabbit pAb
Accession(Primary)	P42224
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; IP-2-5 ug/mg lysate; ELISA-1:10000

Target information (P42224)

Gene Name	STAT1
About Gene Symbol	STAT1 (Signal Transducer And Activator Of Transcription 1) is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell



Western blot analysis of Phospho-STAT1 (Tyr701) in various lysates using Phospho-STAT1 (Tyr701) antibody.