

Phospho-MAP3K8/COT (Thr290) Rabbit pAb


ABZ4412

Product Overview

Name	Phospho-MAP3K8/COT (Thr290) Rabbit pAb
Accession(Primary)	P41279
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (P41279)

Gene Name	MAP3K8
About Gene Symbol	



Western blot analysis of Phospho-COT (Thr290) in 293 lysates treated with UV using Phospho-COT (Thr290) antibody. The lane on the right is blocked with the Phospho- peptide.