

Phospho-STAT3 (Tyr705) (5E9) Mouse mAb


ABZ4413

Product Overview

Name	Phospho-STAT3 (Tyr705) (5E9) Mouse mAb
Accession(Primary)	P40763
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using Phospho-STAT3 (Tyr705) (5E9) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.