

Phospho-STAT3 (Tyr705) Rabbit pAb


ABZ4415

Product Overview

Name	Phospho-STAT3 (Tyr705) Rabbit pAb
Accession(Primary)	P40763
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 88 kDa; Observed MW: 88 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; IP-2-5 ug/mg lysate; ELISA-1:10000

Target information (P40763)

Gene Name	STAT3
About Gene Symbol	STAT3 (Signal Transducer And Activator Of Transcription 3) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-STAT3 (Tyr705)
in various lysates using Phospho-STAT3 (Tyr705)
antibody.