

Phospho-TGF beta Receptor I (Ser165) Rabbit pAb

ABZ4420


Product Overview

Name	Phospho-TGF beta Receptor I (Ser165) Rabbit pAb
Accession(Primary)	P36897
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:2000; IF-1:50-1:300; IHC-1:50-1:100

Target information (P36897)

Gene Name TGFBRI

About Gene Symbol



Western blot analysis of Phospho-TGF beta Receptor I (Ser165) in mouse LIVER lysates using Phospho-TGF beta Receptor I (Ser165) antibody.

