

Phospho-VEGF Receptor 2 (Tyr1175) Rabbit pAb

ABZ4421

Product Overview


Name	Phospho-VEGF Receptor 2 (Tyr1175) Rabbit pAb
Accession(Primary)	P35968
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:10000

Target information (P35968)


Gene Name

KDR

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using Phospho-VEGF Receptor 2 (Tyr1175) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of Phospho-VEGF Receptor 2 (Tyr1175) in MCF-7 cells using Phospho-VEGF Receptor 2 (Tyr1175) antibody, and DAPI (blue).