

Phospho-IRS1 (Ser639) Rabbit pAb


ABZ4424

Product Overview


Name	Phospho-IRS1 (Ser639) Rabbit pAb
Accession(Primary)	P35568
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 132 kDa; Observed MW: 150 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P35568)

Gene Name	IRS1
About Gene Symbol	



Western blot analysis of Phospho-IRS1 (Ser639) in various lysates using Phospho-IRS1 (Ser639) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma, using IRS1 (Phospho-Ser639) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.