

Phospho-beta Catenin (Ser37) Rabbit pAb

ABZ4425


Product Overview

Name	Phospho-beta Catenin (Ser37) Rabbit pAb
Accession(Primary)	P35222
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (P35222)

Gene Name	CTNNB1
About Gene Symbol	

(kD)



Western blot analysis of Phospho-beta Catenin (Ser37) in various lysates using Phospho-beta Catenin (Ser37) antibody.