

Phospho-VE Cadherin (Tyr731) Rabbit pAb


ABZ4426

Product Overview

Name	Phospho-VE Cadherin (Tyr731) Rabbit pAb
Accession(Primary)	P33151
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 88 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P33151)

Gene Name	CDH5
About Gene Symbol	



Western blot analysis of Phospho-Cadherin 5 (Tyr731) in HepG2 lysates using Phospho-Cadherin 5 (Tyr731) antibody.