

Phospho-AKT1/2/3 (Thr308) Rabbit pAb


ABZ4430

Product Overview

Name	Phospho-AKT1/2/3 (Thr308) Rabbit pAb
Accession(Primary)	P31749
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Fruitfly
Calculated MW	Calculated MW: 56 kDa; Observed MW: 56 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P31749)

Gene Name	AKT1
About Gene Symbol	AKT1 (AKT Serine/Threonine Kinase 1) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-AKT (Thr308) in COLO205 lysates using Phospho-AKT (Thr308) antibody.