

Phospho-ERK1/2 (Thr202/Thr185) Rabbit mAb


ABZ4434

Product Overview


Name	Phospho-ERK1/2 (Thr202/Thr185) Rabbit mAb
Accession(Primary)	P27361
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44,42 kDa; Observed MW: 44,42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (P27361)

Gene Name	MAPK1/MAPK3
About Gene Symbol	MAPK3 (Mitogen-Activated Protein Kinase 3) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-ERK1 (Thr202)/ERK2 (Thr185) in HeLa, HeLa+TPA, C6, C6+TPA lysates using Phospho-ERK1/2 (Thr202/Thr185) antibody.



Western blot analysis of Phospho-p44/42 MAPK (ERK1/2) (Thr202) in C6, HeLa lysates using Phospho-p44/42 MAPK (ERK1/2) (Thr202) antibody.