

Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) Rabbit


ABZ4435

Product Overview

Name	Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) Rabbit mAb
Accession(Primary)	P27361
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 44,42 kDa; Observed MW: 44,42 kDa
Recommended Dilution	WB-1:500-1:50000; IP-1:50-1:100

Target information (P27361)

Gene Name	MAPK1/MAPK3
About Gene Symbol	MAPK3 (Mitogen-Activated Protein Kinase 3) is a mitogen-activated protein kinase involved in MAPK/ERK signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Erk1 (pT202/pY204)/Erk2 (pT185/pY187) in C6, HeLa lysates using Phospho-ERK1/2 (Thr202/Tyr204)/(Thr185/Tyr187) antibody.