

# Phospho-Phospholamban (Ser16) Rabbit mAb

ABZ4436


## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-Phospholamban (Ser16) Rabbit mAb     |
| <b>Accession(Primary)</b>   | P26678                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse                                 |
| <b>Calculated MW</b>        | Calculated MW: 6 kDa; Observed MW: 12,24 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                             |

## Target information (P26678)

**Gene Name**

PLN

**About Gene Symbol**

Western blot analysis of Phospho-lamban (Phospho-Ser16) in mouse heart lysates using Phospho-Phospho-lamban (Ser16) antibody.

