

Phospho-Glucocorticoid Receptor (Ser226) Rabbit mAb


ABZ4447

Product Overview


Name	Phospho-Glucocorticoid Receptor (Ser226) Rabbit mAb
Accession(Primary)	P04150
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 86 kDa; Observed MW: 94 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000

Target information (P04150)

Gene Name	NR3C1
About Gene Symbol	NR3C1 (Nuclear Receptor Subfamily 3 Group C Member 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms.



Western blot analysis of Phospho-Glucocorticoid Receptor (Ser226) in HeLa, A549, C6 lysates using Phospho-Glucocorticoid Receptor (Ser226) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Glucocorticoid Receptor (Phospho-Ser226) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.